



MACHEREY-NAGEL
Filtration · Rapid Tests · Water Analysis · Chromatography · Bioanalysis
Filtration · Schnellteste · Wasseranalytik · Chromatographie · Bioanalytik

[Editors](#) [English / Deutsch](#)

[Home](#) [SPE](#) [HPLC](#) [TLC](#) [GC](#) [Flash/Adsorbents](#) [Vials/Closures](#) [Filters](#) [Service](#) [Applications](#) [WebShop](#)

[back to summary of RP phases](#)

CHROMABOND® C₁₈ Hydra

Octadecyl modified silica phase for SPE of polar analytes

Features

- Base material silica, pore size 60 Å, particle size 45 µm, specific surface 500 m²/g, pH stability 2–8
- Special octadecyl phase for polar analytes, not endcapped, carbon content 15 %
- Recommended application:
more polar compounds like pesticides and their polar degradation products, phenols, phenoxycarboxylic acids, nitroaromatics, pharmaceuticals

For applications click here.

Ordering information

	Volume	Adsorbent weight				Pack of		
CHROMABOND® C₁₈ Hydra polypropylene columns								
		50 mg	100 mg	200 mg	500 mg	1 g	2 g	3 g
	1 mL	730294	730295					100
	3 mL			730296	730297	730298		50
	6 mL				730299	730300	730301 730302	30
	CHROMABOND® C₁₈ Hydra glass columns							
				200 mg	500 mg	1 g		
	3 mL		730296G	730297G	730298G			50
	6 mL			730299G	730300G			30
CHROMABOND® LV-C₁₈ Hydra								
				200 mg				
	15 mL			732295				30
CHROMAFIX® C₁₈ Hydra cartridges								
	Size / Adsorbent weight (average)	S	M	L				
		270 mg	530 mg	950 mg				
		731730	731731	731732				50
CHROMABOND® MULTI 96 C₁₈ Hydra								
						96 x 100 mg		
						738294.100M		1
CHROMABOND® C₁₈ Hydra adsorbent								
						730628		100 g

Copyright 2015 by MACHEREY-NAGEL GmbH & Co. KG Imprint / Impressum Privacy Statement Terms Of Use